

Embryo shipping plate

Printing stringency: Low

Our condition

- Layer height: 0.1
- Infill: 20%
- Support: No
- 200 °C, with PLA

Features

- Single print
- 100 wells (10 x 10)
- Embryos are safely preserved with moisture for at least an overnight shipping, which can be extended to >3 days if seal the cap with Scotch tape. Heat pad & plastic bag optional.
- Four components:
 - **Embryo_plate_grid** (need to use transparent filament for making embryos in well visible)
 - **Embryo_plate_cap**
 - **Embryo_plate_band1**
 - **Embryo_plate_band2**

v1.0: Bod Owens, 2021/11/01



